

CO adsorption on HKUST-1 at RT in drifts

Date: 2023-09-01

Tags: HKUST-1 DRIFTS

Created by: Shuang Chen

Date: 28.07.2023

Sample: HKUST-1 (sample Nr. 22 in the Database)

General description of the experiment:

To evaluate the adsorption bands of CO (1 atm) on HKUST-1 at RT (30 °C)

Detailed procedure:

Preparation: Put the water trap and trap containing active carbon into a oven (425 °C) overnight to remove all water and carbonyls.

----- **KBr powders,**

Run a background in an Ar (1 atm, 50 ml/min) atmosphere

Background parameter: 256 scans, 2 cm⁻¹

heat 500 °C in Ar for 1 hour and cooled down in Ar to 30 °C

Run a background in an Ar atmosphere to compare the water/solvent content.

At rt switch to CO (1 atm, 200 ml/min), measure CO adsorption as a function of time up to 1 h.

parameter: 8 scans, 2cm⁻¹

Switch to Ar, flush out CO.

----- **Change sample to HKUST-1,**

Run a background in an Ar atmosphere (1 atm, 50 ml/min)

parameter: 256 scans, 2cm-1

heat 100 °C in Ar for 1 hour and cooled down in Ar to 30 °C

Run a background in an Ar atmosphere to compare cleanliness of the sample

At rt switch to CO (1 atm, 200 ml/min), measure CO adsorption as a function of time up to 1 h.

parameter: 8 scans, 2cm-1

Switch to Ar, flush out CO.

Clean the sample chamber.

Settings

File settings:

Computer: for drift experiment

Folder: D:\data\shuangchen\2023---HKUST-1

Spectra names: SC0001---

Opus settings:

Resolution	2
Phase resolution	4
Acquisition Mode	Double Sided, Forward-Backward
Aperture	8mm
Mode	Absorbance
Detector	LN MCT Narrow [Internal Pos.2]

Gain	Automatic
ADC-counts	
Polarization	no
Scanner velocity	40

Protocol of the experiment:

Action/Spectrum	Name	Temperature °C	Time	Comments
Action	Take traps out of the oven and cool down to RT			
Action	Load the KBr into the DRIFTS cell			
Action	Check for leaks using Ar (1 atm/ 50 ml/min)			
Action	test an initial background of KBr in Ar			
Action	heat the sample in Ar at 500 °C for 1h (to clean sample)			
Action	Cool down the sample in Ar to 23 °C			
Action	test a background of KBr in Ar (after heating to check cleanliness)			
Action	Run the experiment of CO (1 atm/ 200 ml/min) adsorption for 1 hours at 23 °C			
Action	Switch to Ar flush 10min to remove CO			
Action	Change sample to HKUST-1 and put it into the DRIFTS cell			
Action	test a background of HKUST-1 in Ar (same parameter)			
Action	Heat the HKUST-1 in Ar at 100 °C for 1h			
Action	Cool down the HKUST-1 in Ar to 30 °C			
Action	test a background of HKUST-1 in Ar after heating to check cleanliness			
Action	Switch to CO and run the experiment of CO adsorption on HKUST-1 for 1 hours at 30 °C			
Action	Flow Ar for 10min to remove CO			
Action	Remove sample and clean out DRIFTS cell with water			

Results :

Comments:

Steps

Put traps into the oven / clean at 425°C over night

Load the KBr into the DRIFTS cell

Check for leaks using Ar (1 atm/ 50 ml/min)

test an initial background of KBr in Ar

Heat the sample in Ar at 500 °C for 1h (to clean sample)

Cool down the sample in Ar to 30 °C

test a background of KBr in Ar (after heating to check cleanliness)

Run the experiment of CO (1 atm/ 200 ml/min) adsorption for 1 hours at 30 °C

Switch to Ar flush 10min to remove CO

Change sample to HKUST-1 and put it into the DRIFTS cell

test a background of HKUST-1 in Ar (same parameter)

Heat the HKUST-1 in Ar at 100 °C for 1h

Cool down the HKUST-1 in Ar to 30 °C

test a background of HKUST-1 in Ar after heating to check cleanliness

Switch to CO and run the experiment of CO adsorption on HKUST-1 for 1 hours at 30 °C

Flow Ar for 10min to remove CO

Remove sample and clean out DRIFTS cell with water



Unique eLabID: 20230804-648209758ead69d055605cf8965d2bd8493d2230
Link: <https://ifgselabftw.ifg.kit.edu/experiments.php?mode=view&id=2417>