

Treating Ceria Powder at High Temperatures in an Oxygen Atmosphere

Date: 2022-12-21

Tags: DRIFTS CO CeO₂ powder High temperature

Created by: Lachlan Caulfield

Goal :

To Investigate the cleaning effect of O₂ at elevated temperatures.

Procedure :

Background run in oxygen before heating (2 cm⁻¹) - LC0028.0

Heat sample to 1000 °C for 1 hour

Take second background after cooling

Increase temperature to 550 °C and flow CO 200 ml/min

Use background (2 cm⁻¹) - LC0028.1

Heating 1 hour in O₂ (start 09:45)

Cool down in O₂ (start 11:45)

Flow in CO at 75 °C, after O₂ treatment - LC0028.0000-0139

heat to 100 °C at LC0028.0140

heat up carbonyl trap to 400 °C using the heat gun

increase water flow into the cell

LC0028.0141-0173 (14:38) gun switch off

Results :

It was found that the carbonyl trap stores water over time, therefore it must be replaced after each experiment, otherwise the water content will be too high for CO adsorption on metal oxides.

(Data saved - DRIFTS PC; Folder - Data --> L Caulfield; File name - 20221221_CeO₂_O₂_treatment_xTemp)



Unique eLabID: 20230704-676b28977ab141c51109a66efabca087be74d735
Link: <https://ifgselabftw.ifg.kit.edu/experiments.php?mode=view&id=2265>