

Ceria Cleaned with O₂ at High Temperatures and Exposed to CO at 1 Atm

Date: 2023-02-15

Tags: DRIFTS CO CeO₂ powder High temperature Cleaned

Created by: Lachlan Caulfield

Goal :

To repeat earlier experiments that were cleaned in O₂ prior to CO exposure, however, this time the pure CO is used.

Procedure :

New sample, heated overnight 225 °C

LC0040.0 - in Ar atmosphere

Flush for 5 additional minutes to remove residual CO gas

LC0040.0000-0002

LC0040.1 - in O₂ atmosphere

Heat sample to 1000 °C, 2 hours, start time 09:36, reached temp (764 °C) at 09:45

Cool down start at 11:45

Final temp 24 °C

Second background used --> Use background from before heating to stay consistent with 20230215

Switch off oxygen and flow through CO - LC0040.0003-502

LC0040.0503-612, 256 scans

water trap warm, change to Argon, 50 ml/min

LC0040.0613 at 256 scans

Results :

Data collected was the same as collected earlier when the ceria sample is cleaned in O₂, the experimental procedure is reproducible and now a subtraction method is required.

(Data saved - DRIFTS PC; Folder - Data --> L Caulfield; File name - 20230215_CeO₂-HightT_O₂_pressure_1bar)



Unique eLabID: 20230704-d91787b5d698d920488621646365d1bf8f7b079f

Link: <https://ifgselabftw.ifg.kit.edu/experiments.php?mode=view&id=2269>