

# Reduction of Ceria after 1 Hour of CO Exposure

**Date:** 2022-11-30

**Tags:** CeO<sub>2</sub> DRIFTS CO room temperature

**Created by:** Lachlan Caulfield

---

## Goal :

To understand the extent of reduction over a period of 1 hour and compare it to 2 hours of CO exposure.

## Procedure :

Run a background in an O<sub>2</sub> atmosphere preheating with a resolution of 2 cm<sup>-1</sup> LC0011.0  
Heating in O<sub>2</sub> for 1 hour 1000 °C

Initial sample is lost - Try again (Blown by gas)

Run second preheating background also in an O<sub>2</sub> atmosphere with a resolution of 2 cm<sup>-1</sup> LC0011.1

-----

Heat for 78 minutes in O<sub>2</sub>, at 1000 °C

Cool the sample down to room temperature

T < 30 °C in O<sub>2</sub> batch, change over time - LC0011.0000-0001

Go to bypass and run CO through bypass during cool down phase (~ 30-60 min) start 13:55 - 14:32 (37)  
When at room temperature switch CO to sample from the bypass (26 °C) - LC0011.0002-0012

Run the experiment for 1 hour at 100 °C (ref to 13 change during heating) - LC0011.0013-0047

Heater turned off 15:56

LC0011.0048-0059 - cool down, measure at rt, probably hydroxylation (ref to 48, look to cool down)

(start 24 °C) (ref to 12 ... before after heating - LC0011.0060-0074

no real change anymore

flush with O<sub>2</sub>, change setup, heat to 200 °C, change to CO, with LN<sub>2</sub> trap,

Sample lost, try again

---

LC0011.0079-0081 - mit LN<sub>2</sub> trap, 200 °C

switch to CO

LC0011.0082-0084 - at 200 °C

cool down in CO, with LN<sub>2</sub> trap

## Results :

No conclusive signs of reduction, the sample may not have been properly cleaned, LN<sub>2</sub> trap is tested to see if it can reduce the water content. Water leads to poor adsorption.

**(Data saved - DRIFTS PC; Folder - Data --> L Caulfield; File name - 20221130\_CeO2\_60\_red)**



Unique eLabID: 20230703-ed192011c4baf9fc1eae3d19049a0e063f80b991  
Link: <https://ifgselabftw.ifg.kit.edu/experiments.php?mode=view&id=2258>