

# Reduction of Ceria at Room Temperature and 150 °C in a CO Atmosphere

**Date:** 2023-04-14

**Tags:** DRIFTS CO room temperature CeO<sub>2</sub> powder High temperature

**Created by:** Lachlan Caulfield

---

## Goal :

To see the effect CO has on the bands when going from room temperature to 150 °C in a CO atmosphere.

## Procedure :

Background run in O<sub>2</sub> at room temperature before heating - LC0064.0

Remove excess CO  
LC0064.1-2

-----

Heat sample in O<sub>2</sub> for 2 hours 550 °C (Start 09:30)  
Cool down to 150 °C

-----

Run background at 150 °C - LC0064.3  
Introduce CO - In flow  
LC0064.0000-0705

-----

Change samples

-----

Before heating in O<sub>2</sub> background - LC0065.0  
Heat for 2 hours in O<sub>2</sub>  
Cool to rt before CO  
Measure a background after heating in O<sub>2</sub> - LC0065.4  
LC0065.0000-0300

switch to 2ml/min CO

## Results :

Data was collected for both room temperature and 150 °C and can be used for subtraction.

(Data saved - DRIFTS PC; Folder - Data --> L Caulfield; File name - 20231304\_CeO2\_rt\_150)



Unique eLabID: 20230706-e27d5975efc1497393ce043a2805ec29c3ec4ed5

Link: <https://ifgselabftw.ifg.kit.edu/experiments.php?mode=view&id=2305>