

# Reduction of Ceria after 2 Hours of CO Exposure

**Date:** 2022-11-29

**Tags:** CeO2 DRIFTS CO room temperature CeO2 powder

**Created by:** Lachlan Caulfield

---

## Goal :

To observe reduction over a long period of time on a cleaned ceria powder.

## Procedure :

Treated 1 hour in oxygen at 1000 °C

background run in an O<sub>2</sub> atmosphere with 2 cm<sup>-1</sup> resolution

-----  
Highly impure reading leaving to run over time to see if it disappears  
CO 1 bar 200 ml/min

Removal of bulk impurities - LC0010.0000-0027

Room temperature in CO - LC0010.0028-0037

100 CO 120 min (start 13:37 end 15:37) - LC0010.0038-0103

Cool down in CO (end 27 °C) - LC0010.0104-0113

Bands at ~2175 cm<sup>-1</sup> and 2168 cm<sup>-1</sup>, lose over time with increase of 3695 cm<sup>-1</sup> -> switch CO for OH, also small changes in carbonates

Room temperature, CO, hydroxylation? 30 min, then transient DRIFTS - LC0010.0113-0151

change to 50 ml/min Ar, measure 2 cm<sup>-1</sup>, 24 scans (6 s > spectrum recording time > 5 s), cool before (22 °C)

50 ml/min Ar, over night LC0010.0152-1151

Recorded a spectra the next day to see the change over time - LC0010.0152

## Results :

Bands were observed under the CO gas phase at 2175 cm<sup>-1</sup> and 2168 cm<sup>-1</sup> which are related to CO ads on ceria. Interesting side reactions occur with OH and carbonate formations.



Unique eLabID: 20230703-1c96507309505389dc4d1d4fb38e64a985847076

Link: <https://ifgselabftw.ifg.kit.edu/experiments.php?mode=view&id=2257>